
One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) For self-test

Product Pictures

➤ Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (For Self-test)


1 Test /Kit

5 Tests /Kit

User Manual

Biohazard Sample Bag

Sterile Swab

SARS-CoV-2 Antigen Test Card

Extraction Tube with Sample Extraction Solution and Tip

Specifications

➤ Info. of the Test Kit and Export Packing Cartons

Product Name	Specifications	Size of Kit (L x W x H)	Weight/Kit	Kit Quantity/Carton	Size of Carton (L x W x H)	Weight per Carton
One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)	1 Test/Box	7cm x 1.8cm x 13 cm	25.5 g	300 Boxes/Carton	49.5cm x 40.5cm x 30 cm	8.6 kg/Carton
	5 Tests/Box	7cm x 5.3cm x 13 cm	76 g	120 Boxes/Carton	55.5cm x 41cm x 30 cm	10.44 kg/Carton

Export Packing Cartons

Packing carton for 5 T/kit
Size: 55.5*41*30 cm

Packing carton for 1 T/kit
Size: 49.5*40.5*30 cm

Brochure


STOCK CODE 603387

One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) • For self-test

CE Marked


*Rapid, Convenient, Easy and Reliable
Detection of COVID-19*


Intended Use

One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) is intended for the qualitative detection of SARS-CoV-2 antigens in human nasal swab samples.

This test is suitable for medical laypersons as a self-test at home or at work.

Product Components


SARS-CoV-2 Antigen Test Card

Biohazard Sample Bag

Extraction Tube with Sample Extraction Solution and Tip


Sterile Swab

User Manual

COVID-19 Diagnostic Testing

- **PCR-based tests** can detect small amounts of viral genetic materials.
- **Antigen tests** detect the presence of viral proteins and can return positive results when a person is most infectious.
- - - **Antibody tests** detect the body's immune response to the virus.

Features

Non-invasive sampling
(Sample type: nasal swab)

Read test results visually.
Do not require test equipment.

Early detection of SARS-CoV-2 infection

Rapid test.
Test result available in 10-15 min.

Simple operation,
easy to learn and use

Room temperature storage (4-30°C)

When to Use the Test Kit?

Use this test:

- ✔ If you want to test yourself.
- ✔ If you have symptoms similar to COVID-19, such as: E.g. headache, fever, cough, sore throat, loss of sense of smell or taste, shortness of breath, muscle pain.
- ✔ If you are concerned about whether you are infected with COVID-19.
- ✔ Use of the test by persons under 18 years of age only under the supervision of an adult.

Please Scan the QR Code to Get
Operation Video and Multi-language User Manual

Operation
Video


Multi-language
User Manual

Test Procedure

Nasal Swab
Sampling

>>

Insert the collection tip
of swab into one nostril.


Rotate ≥ 4 times
(15 s in total)


Repeat with the
same swab.

Sample Pretreatment >>

Test >>

Test Results

Positive

Positive

Negative

Invalid 1

Invalid 2

Specifications

Product Name	One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)
Test Item	SARS-CoV-2 Antigen
Package	1 T/kit, 5 T/kit, 7 T/kit, 25 T/kit
Product Code	CG20615/ CG206155/ CG206157/ CG2061525
Test Time	10-15 min
Storage Condition	4-30°C
Shelf Life	24 months
Recommended Test Temperature	23-25°C

Application Scenarios

Home

School

Work Place

Dormitory

Nursing Home

Cruise ship

Airport

Theater

Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China

Tel: +86-25-68568508/68568594

Fax: +86-25-68568500

E-mail: sales@getein.com.cn; overseas@getein.com.cn

Web: en.bio-gp.com.cn

CMC Medical Devices & Drugs S.L.

Add: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain

Tel: +34951214054

User Manual

➤ Please scan the QR code to get the user manual of your language.

Languages:

- English
- French
- German
- Danish
- Dutch
- Spanish
- Portuguese
- Italian
- Norwegian
- Finnish
- Swedish


Multi-language user manual

One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)

User Manual for self-testing

INTENDED USE
One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) is intended for the qualitative detection of SARS-CoV-2 antigens in human nasal swab samples. This test is used for individuals suspected of COVID-19 within the first seven days of the onset of symptoms, such as headache, fever, cough, sore throat, loss of the sense of smell or taste, shortness of breath, muscle pain. Meanwhile the test can also be applied for individuals without symptoms. Results are for the identification of SARS-CoV-2 antigen. Positive results indicate the presence of SARS-CoV-2 antigens, but individual history and other diagnostic information is necessary for determine infection status. Negative results do not rule out SARS-CoV-2 infection. Negative results for individuals with symptoms similar to COVID-19 infection for more than seven days should be treated as negative possibly. If necessary, it should be confirmed by molecular assay. One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) is intended to be used to help the diagnosis of SARS-CoV-2 infection. This test is used for self-testing.


PREPARING THE TEST


- Check integrity of the out package, components and the expiration date.
- Read the user manual before starting the test. Check operation video for more help.
- Open the pouch, check the result window and sample well (S).

TEST PROCEDURE


- Remove the kit from the top of the extraction tube with sample extraction solution.

TEST RESULTS

SPECIMEN COLLECTION

STORAGE AND STABILITY
Store the test kit at 4-30°C with a valid period of 24 months. Use the test card within 1 hour once the foil pouch is opened.

PRINCIPLE
The test uses anti-SARS-CoV-2 nucleocapsid protein (N protein) monoclonal antibody conjugated with colloidal gold coated on the sample pad, and another anti-SARS-CoV-2 N protein monoclonal antibody (I) coated on test line. After the samples have been applied to the test strip, the colloidal gold-labeled anti-SARS-CoV-2 N protein monoclonal antibody (I) bind with SARS-CoV-2 antigens in sample and form marked antigen-antibody complexes. These complexes move to the test card detection zone by capillary action. Then marked antigen-antibody complexes will be captured on test line by anti-SARS-CoV-2 N protein monoclonal antibody (I). The color intensity of each test line increases in proportion to the amount of SARS-CoV-2 antigen in sample.

PRECAUTIONS
1. Always keep out of the reach of children. Small parts of the kit can be a choking hazard.
2. Sample extraction solution is a phosphate buffer contained low concentration of sodium chloride, Tween, hexadecyl trimethyl ammonium bromide and sodium azide. If extraction solution splashes your body or into eyes, please wash with water.

LIMITATIONS
1. False-negative result may occur if the level of antigen in sample is below the detection limit of the test or the sample was collected incorrectly.
2. Clinical diagnosis and treatment cannot be made without consulting with the physician.
3. A negative result from an individual have symptoms similar to COVID-19 beyond seven days should be treated as negative possibly. If necessary, confirmed with the molecular assay.
4. The product One Step Test For SARS-CoV-2 Antigen (Colloidal Gold) showed no drop off in sensitivity when compared with the wild type with respect to the following variants: VOC1 UK, Alpha, VOC2 South Africa, Beta, VOC3 Brazil, Gamma, VO1 America Iota and VO2 India Kappa. We will keep evaluating the spread of new variants.

PERFORMANCE CHARACTERISTICS

1 Limit of Detection (LoD)
The LoD for nasal swab was established using heat-inactivated SARS-CoV-2 isolate strain. The strain was spiked with negative human nasal swab into a series of concentrations. The estimated LoD found from the initial twofold serial dilution test was confirmed by testing 20 replicates. The confirmed LoD for nasal swab was 200 TCID₅₀/mL.

2 Clinical Agreement Study
The clinical performance of One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) was evaluated by testing a total of 480 nasal swab samples. It is compared to the results of RT-PCR assays. Overall study results were shown in the tables below.

GeneX kit	BGI RT-PCR kit		
	Positive	Negative	Subtotal
Positive	165	4	169
Negative	5	396	401
Subtotal	170	400	480

Positive percent agreement (Diagnostic sensitivity) = 165 / (165 + 5) × 100% = 97.06% (95% CI: 93.30%-99.74%)

Negative percent agreement (Diagnostic specificity) = 306 / (306 + 4) × 100% = 98.71% (95% CI: 96.27%-99.59%)
Total percent agreement = (165 + 306) / 480 × 100% = 98.13% (95% CI: 96.48%-99.11%)

3 Analytical Specificity

3.1 Cross-Reactivity & Microbial Interference
Each organism and virus was tested in triplicate in the absence and presence of SARS-CoV-2 respectively. According to the test results, there was no cross-reactivity with the following viruses or organisms.

Viruses or organisms	Concentration
Human coronavirus 229E	1 x 10 ⁷ PFU/mL
Human coronavirus OC43	1 x 10 ⁷ PFU/mL
Human coronavirus NL63	9.87 x 10 ⁷ PFU/mL
MERS coronavirus	7930 PFU/mL
Adenovirus (e.g. C1 Ad. 71)	1 x 10 ⁷ PFU/mL
Human Metapneumovirus (hMPV)	1 x 10 ⁷ PFU/mL
Parainfluenza virus Type 1	1 x 10 ⁷ PFU/mL
Parainfluenza virus Type 2	1 x 10 ⁷ PFU/mL
Parainfluenza virus Type 3	1 x 10 ⁷ PFU/mL
Parainfluenza virus Type 4a	1 x 10 ⁷ PFU/mL
Influenza A	1 x 10 ⁷ PFU/mL
Influenza B	2.92 x 10 ⁷ PFU/mL
Enterovirus	1 x 10 ⁷ PFU/mL
Respiratory syncytial virus	1 x 10 ⁷ PFU/mL
Rhinovirus	4.17 x 10 ⁷ PFU/mL
Haemophilus influenzae	1 x 10 ⁷ CFU/mL
Streptococcus pneumoniae	1 x 10 ⁷ CFU/mL
Streptococcus pyogenes	1 x 10 ⁷ CFU/mL
Cardiobacterium	1 x 10 ⁷ CFU/mL
Pooled human nasal wash	14% v/v
Bordetella pertussis	1 x 10 ⁷ CFU/mL
Mycoplasma pneumoniae	1 x 10 ⁷ CFU/mL
Chlamydia pneumoniae	1 x 10 ⁷ CFU/mL
Legionella pneumophila	1 x 10 ⁷ CFU/mL
Mycobacterium tuberculosis	1 x 10 ⁷ CFU/mL
Pneumocystis jirovecii	1 x 10 ⁷ CFU/mL
Pseudomonas Aeruginosa	1 x 10 ⁷ CFU/mL
Staphylococcus Epidermidis	1 x 10 ⁷ CFU/mL
Streptococcus Salivarius	1 x 10 ⁷ CFU/mL

3.2 Interferences


The potentially interfering substances that may be found in the upper respiratory tract in symptomatic subjects (including over the counter medications). No false positive or false negative results were seen at the following concentrations.

Potential Interfering Substances	Concentration
Blood (human)	5%
Mucin	5 mg/mL
Nasal mucus (human)	5% v/v
CVS Nasal Drops (phenylephrine)	15% v/v
Afrin (Oxymetazoline)	15% v/v
CVS Nasal Spray (Cromolyn)	15% v/v
Zicam Cold Remedy	5% v/v
Homeopathic (Abulad)	10% v/v
Sore Throat Phenol Spray	15% v/v
Tobramycin	3.3 mg/dL
Mupirocin	0.15 mg/dL
Fluticasone	5% v/v
Tamiflu (Oseltamivir phosphate)	500 mg/dL

Biotin	0.35 mg/dL
Methanol	0.15% v/v
Diphenhydramine	0.071 mg/dL
Dextromethorphan	0.00156 mg/dL
Dexamethasone	1.2 mg/dL

4 Precision
For repeatability study, the agreement percent of both negative samples and positive samples are 100%. For reproducibility study, the agreement percent of both negative samples and positive samples are 100%.

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Certificates and Clinical Evaluation

➤ Clinical Agreement Study

Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (For Self-Test)	
Sensitivity	97.06%
Specificity	98.71%
Total Percent Agreement	98.13%

The clinical performance of One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) was evaluated by testing a total of 480 nasal swab samples. It is compared to the results of RT-PCR assays. Overall study results were shown in the tables below.

Total	BGI's RT-PCR kit		
	Positive	Negative	Subtotal
Getein's kit	165	4	169
	5	306	311
	170	310	480

Positive percent agreement (Diagnostic sensitivity) = $165 / (165 + 5) \times 100\% = 97.06\%$ (95% CI: 93.30%-98.74%)

Negative percent agreement (Diagnostic specificity) = $306 / (306 + 4) \times 100\% = 98.71\%$ (95% CI: 96.73%-99.50%)

Total percent agreement = $(165 + 306) / 480 \times 100\% = 98.13\%$ (95% CI: 96.48%-99.01%)

**CE Certificate of Getein One Step Test
for SARS-CoV-2 Antigen (Colloidal
Gold) (For Self-Test)**

CERTIFICATE

EC Certificate No. 1434-IVDD-447/2021
EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices

Polish Centre for Testing and Certification certifies
that manufactured by:

GETEIN Biotech, Inc.

Nanjing, ul. Bofu Road, Luhe District 9, China

in vitro diagnostic medical devices
for self-testing

One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)

Ref. codes: CG20615, CG206152, CG206153, CG206155, CG206156, CG206157, CG206158, CG206159,
CG2061510, CG2061511, CG2061512, CG2061515, CG2061520, CG2061525

In terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 30.07.2021 to 27.05.2024

The date of issue of the Certificate: 30.07.2021



The date of the first issue of the Certificate: 30.07.2021

Issued under the Contract No. MD-66/2021
Application No: 142/2021
Certificate bears the qualified signature.
Warsaw, 30.07.2021
Module A1

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Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (For Self-Test) has been approved by Malaysia MDA

➤ **Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (For Self-Test)** has been approved by Thailand FDA

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Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) has passed the phase 3a validation of British DHSC

GOV.UK Topics Government activity

→ Coronavirus (COVID-19) | Guidance and support

Home > Coronavirus (COVID-19) > Testing for coronavirus (COVID-19) > Assessment and procurement of coronavirus (COVID-19) tests

Department of Health & Social Care

Guidance

Outcome of the evaluation of rapid diagnostic assays for specific SARS-CoV-2 antigens (lateral flow devices)

Updated 11 October 2021

Background

Since its establishment in August 2020, the joint PHE Porton Down and University of Oxford SARS-CoV-2 lateral flow antigen test validation cell has evaluated over 140 lateral flow devices that have been referred by the Department of Health and Social Care (DHSC).

Approximately 30% of the tests that were referred for validation met the standards for phase 2 validation, which are set out in the [protocol](#). PHE Porton Down subsequently performed phase 3 testing to assess whether the lateral flow devices that passed phase 2 displayed performance characteristics desirable for mass population, community-based testing. The desirable performance characteristics are:

- very high specificity
- very high sensitivity against viral loads associated with infectiousness

The lateral flow devices that display the desirable performance characteristics are summarised in table 1 below.

Getein Biotech One Step Test for SARS-CoV-2 Antigen	Pass	10 March 2021
-----------------------------------------------------	------	---------------

CE Common List

EUROPEAN COMMISSION
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY
Public health, country knowledge, crisis management
Health Security

**EU health preparedness:
A common list of COVID-19 rapid antigen tests,
including those of which their test results are mutually
recognised, and a common standardised set of data to
be included in COVID-19 test result certificates**

Agreed by the Health Security Committee

This document was agreed by the HSC on 17 February 2021

Annex I

Common list of rapid antigen tests

*A first update was agreed by the HSC on 10 May 2021
A second update was agreed by the HSC on 16 June 2021*

Annex II

Common standardised data set of to be included in COVID-19 test result certificates

An update to Annex II was agreed by the HSC on 19 March 2021

Getein Biotech, Inc.	One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)	Yes	97.06% sensitivity 98.71% specificity Nasal swab	DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 90% at <Ct30 and 100% at <Ct25)		DE ⁽²⁾		DE ⁽²⁾		Yes (2183)	No
----------------------	-------------------------------------------------------	-----	--------------------------------------------------------	------------------------------------------------------------------------------------------------------------	--	-------------------	--	-------------------	--	------------	----

Getein COVID-19 Test Kits on the Export Whitelist of China

序号 No.	生产企业 Manufacturer	统一社会信用代码 Certificate for Uniform Social Credit	国外注册认证情况 Certificate	省份 Province of China	产品型号 Product Name
42	基蛋生物科技股份有限公司 Getein Biotech, Inc.	913201007360621166	CE	江苏 Jiangsu	One Step Test for Novel Coronavirus(2019-nCoV) IgG antibody (Colloidal Gold) One Step Test for Novel Coronavirus(2019-nCoV) IgM antibody (Colloidal Gold) One Step Test for Novel Coronavirus(2019-nCoV) IgM/IgG antibody (Colloidal Gold) Novel Coronavirus (2019-nCoV) Real-time RT-PCR Kit Akso SARS-CoV-2 Real-time RT-PCR Kit SARS-CoV-2 Antigen Fast Test Kit (Immunofluorescence Assay) One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) One Step Test for SARS-CoV-2 Total Antibody/Neutralizing Antibody (Colloidal Gold) One Step Test for FluA/FluB/SARS-CoV-2 Antigen (Colloidal Gold) SARS-CoV-2 Total Antibody/Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay) FluA/FluB/SARS-CoV-2 Antigen Fast Test Kit (Immunofluorescence Assay) FluA/FluB/SARS-CoV-2 Real-time RT-PCR Kit SARS-CoV-2/VOC-202012/01 Real-time RT-PCR Kit SARS-CoV-2/501.V2 Real-time RT-PCR Kit Novel Coronavirus(2019-nCoV) IgM/IgG antibody Fast Test Kit (Immunofluorescence Assay) SARS-CoV-2 Neutralizing Antibody Fast Test Kit (Immunofluorescence Assay) One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (Saliva)
623	基蛋生物科技股份有限公司 Getein Biotech, Inc.	913201007360621166	Austria BASG Listing	江苏 Jiangsu	One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (Nasal) Selftest

Certificate of ISO 13485:2016

bsi.	 By Royal Charter	
Certificate of Registration		
QUALITY MANAGEMENT SYSTEM - ISO 13485:2016		
This is to certify that:	Getein Biotech, Inc. No.9 Bofu Road Luhe District Nanjing Jiangsu 211505 China	基蛋生物科技股份有限公司 中国 江苏省 南京市 六合区 沿江工业开发区 博富路9号 邮编：211505
Holds Certificate No: MD 728432		
and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:		
Design & Development, Manufacture and Distribution of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), PCR Assay to diagnose infectious disease. Design & Development, Manufacture and Distribution of Analyzers in use of Chemiluminescence Immunoassay, Biochemistry Assay, Point of Care Assay (including Colloidal Gold Assay, Immunofluorescence Assay, Dry Chemistry Assay), PCR Assay to diagnose infectious disease. 研发、生产和销售化学发光法试剂、生化试剂、即时诊断（包括胶体金法、免疫荧光法、干式化学法）试剂、传染病相关PCR分子诊断试剂。 研发、生产和销售用于化学发光法试剂、生化试剂、即时诊断（包括胶体金法、免疫荧光法、干式化学法）试剂、传染病相关PCR分子诊断试剂配套使用的分析仪。		
		
For and on behalf of BSI:	Gary E Slack, Senior Vice President - Medical Devices	
Original Registration Date: 2020-05-29	Effective Date: 2020-07-26	
Latest Revision Date: 2021-09-23	Expiry Date: 2023-07-25	
Page: 1 of 2		
...making excellence a habit. [®]		
<small>This certificate was issued electronically and retains the property of BSI and is bound by the conditions of contract. An electronic certificate can be authenticated online. Printed copies can be validated at www.bsi.com/online (Charter/Directory) or telephone +44 (0) 800 30000. Information and Contact: BSI, John W. Kayeplein 9, 1000 EP Amsterdam The Netherlands. Tel: +31 (0) 20 5450 789 BSI Group The Netherlands B.V., registered in the Netherlands under number 33264281, at John H. Kayeplein 9, 1056 EP Amsterdam, The Netherlands A Member of the BSI Group of Companies.</small>		

Certificate of MDSAP

bsi.

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

Getein Biotech, Inc.
No.9 Bofu Road
Luhu District
Nanjing
Jiangsu
211505
China

Facility ID Number: F004902

Holds Certificate No:

MDSAP 728434

Statement of Conformity: The company listed on this certificate has been audited to and found to conform with the following criteria: ISO 13485:2016 and Australia - Therapeutic Goods (Medical Devices) Regulations, 2002, Schedule 3 Part 1 (excluding Part 1.6) - Full Quality Assurance Procedure; Brazil - RDC ANVISA n. 16/2013, RDC ANVISA n. 23/2012, RDC ANVISA n. 67/2009; Canada - Medical Devices Regulations - Part 1 - SCR 98/282; Japan - MHLW Ministerial Ordinance 169, Article 4 to Article 68, PMD Act; USA - 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 - Subparts A to D

Please see scope page.

For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2020-10-18

Effective Date: 2020-10-18

Expiry Date: 2023-07-25

Page: 1 of 2

MDSAP

BSI Group America Inc. is an MDSAP authorized auditing organization

...making excellence a habit™

Thanks for Your Attention